



## Left-Atrial Appendage Closure Device

- *ultraseal* LAA conforms naturally to any Left Atrial Appendage
- 12 stabilizing anchors
- Fully retrievable and repositionable within the Appendage
- Atraumatic soft lobe allows safe advancement into the LAA



The *ultraseal* LAA device has been developed for the transvenous closure of the Left Atrial Appendage.


The patented double articulating center of the *ultraseal* allows for optimal positioning and repositioning of the device within the Appendage.



## ultraseal LAA Device

Order Number	Lobe/Ball Diameter	Disc Diameter	Delivery Sheath
LAA-16	16mm	22mm	10F
LAA-18	18mm	24mm	10F
LAA-20	20mm	26mm	11F
LAA-22	22mm	28mm	11F
LAA-24	24mm	30mm	11F
LAA-26	26mm	32mm	11F
LAA-28	28mm	34mm	12F
LAA-30	30mm	36mm	12F
LAA-32	32mm	38mm	12F

45° x 45° double curve tip on the Cardia Delivery Sheath allows for easier access to the LAA



## cardia LAA Delivery System

Product Code	Description	Sheath
AC-10-0802	10F Standard Kit	80 cm length with LAA curve
AC-11-0802	11F Standard Kit	80 cm length with LAA curve
AC-12-0802	12F Standard Kit	80 cm length with LAA curve

cardia, Inc. • 2900 Lone Oak Parkway • Eagan, MN 55121 USA • Phone +1-651-691-4100

• Website: [www.cardia.com](http://www.cardia.com)

• LinkedIn: [linkedin.com/company/cardia-inc.](https://www.linkedin.com/company/cardia-inc.)

• Fax +1-651-994-9500

• Twitter: @CardiaInc

mkt.br.2016.02.01.3  
for U.S. export only